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(en)

Biodentine™ is a bioactive dentine substitute from the “Active Biosilicate Technology™” innovation.

#### CONTENT

Assembly kit containing: <ul style="list-style-type: none"> <li>• 5 capsules of 700 mg</li> <li>• 5 single-dose containers of 0.20 ml</li> <li>• 5 spatulas</li> <li>• 1 capsule holder.</li> </ul>	Assembly kit containing: <ul style="list-style-type: none"> <li>• 15 capsules of 700 mg</li> <li>• 15 single-dose containers of 0.20 ml</li> <li>• 15 spatulas</li> <li>• 1 capsule holder.</li> </ul>
---	--

#### COMPOSITION

Powder is composed of tricalcium silicate, zirconium oxide, calcium oxide, calcium carbonate and iron oxides.

Aqueous solution is composed of calcium chloride, polycarboxylate and purified water.

#### INTENDED USE

Bioactive dentine substitute

#### TARGET POPULATION

Biodentine™ is intended to be applied on temporary teeth, on permanent teeth (mature or immature) or on non-vital immature permanent teeth depending on indication. Please refer to the dedicated section (§ Indications for use).

#### INTENDED USER

Dental healthcare professional only

#### INDICATIONS FOR USE

##### In the crown:

- For temporary teeth and permanent teeth (immature or mature)
- Permanent dentine restoration under composites or Inlay/Onlay
- Temporary dentine-enamel restoration

- Restoration of deep and/or large coronal carious lesions (sandwich technique).
- Restoration of cervical or radicular lesions

##### On the pulp:

- For temporary teeth and permanent teeth (immature or mature)
- Pulp capping (direct and indirect)
- Pulpotomy for diagnosed symptoms of reversible pulpitis and irreversibly pulpitis where bleeding is controlled within 5 minutes

##### In the root:

For permanent teeth (immature or mature)

- Repair of root perforations
- Repair of furcation perforations
- Repair of perforating internal resorptions
- Repair of external resorption
- Root-end filling in endodontic surgery (retrograde filling)

For immature permanent teeth:

- Apexification
- Revitalization procedure by means of revascularization of permanent immature teeth with necrotic pulp

#### CONTRAINDICATIONS

- Hypersensitivity to one of the ingredients
- Restoration of large loss of tooth substance and subjected to high stresses
- Pulp space that requires post or core for the final restoration
- Biodentine™ is not to be used for permanent obturation or sealing of root canal.

#### SIDE EFFECTS

Possible hypersensitivity to Biodentine™

#### PROPERTIES

Biodentine™ is a bioactive dentine substitute from the « Active Biosilicate Technology™ » innovation.

- Biodentine™ has mechanical properties similar to the sound dentine and can replace it both in the crown and in the root, without any preliminary conditioning of mineral tissues
- Biodentine™ contains high-purity, monomer-free mineral ingredients and is highly biocompatible.
- Biodentine™ creates optimal conditions for maintenance of pulp vitality, by providing a very tight seal on the dentine surface. It therefore reduces risk of postoperative sensitivity and longevity of restorations in vital teeth.
- Biodentine™ being bioactive, creates optimal environment to build reactionary dentine. Dentine bridges are created faster and are thicker than with similar dental materials and represent necessary condition for optimal pulp healing.
- Biodentine™ offers an initial setting time of 12 minutes from start of mix, for optimal use in the crown.

#### CLINICAL PERFORMANCE

- Preservation of the pulp
- Preservation of the amelo-dentinal structure
- Preservation of the root canal structure
- Root edification
- Sealing ability necessary to create a favorable environment required in revitalization procedure for pulp stem cells

#### CLINICAL BENEFITS

- Preservation of the tooth vitality (except for apexification, perforations and furcations)
- Preservation of the tooth functionality
- Revitalization of a necrotic tooth

#### WARNINGS AND PRECAUTIONS OF USE

Warnings:

- Do not manipulate the device without appropriate personal protective equipment because it may cause eye irritation.
- Use a rubber dam to ensure that the cement is not swallowed.
- Wait until the end of the setting time before performing the permanent enamel restoration or removing the matrix (when used).
- Single-use product: do not use the mixed paste on several patients to avoid cross contamination.
- Biodentine™ powder and Biodentine™ liquid are only intended to be mixed together. Otherwise, product performances are not guaranteed.
- Use the recommended ratio of Biodentine™ powder and Biodentine™ liquid. Otherwise, product performances are not guaranteed and may lead to the failure of the procedure.
- For revitalization procedure, wait until the formation of blood clot before applying Biodentine™.

Caution:

- Water contamination slows the setting of the material. Prevent exposure to water and fluids during the initial setting stage.

Note:

- For revitalization procedure, during the second appointment, it is recommended to use a local anesthetic without vasoconstrictor in order to create a blood clot.

#### USE INSTRUCTIONS (FOR EACH INDICATION)

##### Biodentine™ mixing instructions

- Open a capsule and place it on the white capsule holder.
- Detach a single-dose container of liquid.
- Twist cap to open. Be careful that no drop of liquid falls out of the single dose container.
- Pour 5 drops from the single-dose container into the capsule.
- Close the capsule. Place the capsule on a mixing device, such as Technomix, SYG-200, Tac 400 (Lineatac), Silamat, Cap-Mix, Rotomix, Ultramat etc., at a speed of 4000 – 4200 rotations/min.
- Mix for 30 seconds.

- Open the capsule and check the material's consistency.
- If a thicker consistency is preferred, wait for 30 sec to 1 min before checking again. If after trituration the material has a granular and not a creamy appearance, add an additional drop of liquid and mix in the amalgamator for an additional 10 seconds. In this case, the setting time may be increased.
- Collect Biodentine™ with the spatula supplied in the box. Depending on the desired application, you may handle Biodentine™ with an amalgam carrier, a spatula or a Root Canal Messing Gun.

Rapidly rinse and clean the instruments to remove any residual material.

##### IMMEDIATE ENAMEL RESTORATION:

Assess pulp vitality by the usual tests.

- Isolate the tooth with a rubber dam.
- Remove the infected dentine with a round bur and/or a hand excavator. Leave the affected dentine.
- Adapt a matrix around the tooth if a wall is missing.
- Prepare Biodentine™ as indicated above (Biodentine™ mixing instructions)
- Insert Biodentine™ in the cavity, so that the volume of missing dentine is replaced by the same volume of Biodentine™ avoiding to trap air bubbles. Flatten the material without excessive pressure and ensure good adaptation to the cavity walls and margins.
- Wait until the end of the setting time before performing the permanent enamel restoration. Biodentine™ is compatible with all direct crown restoration techniques and particularly with all types of bonding systems.

##### NON-IMMEDIATE ENAMEL RESTORATION

Assess pulp vitality by the usual tests.

- Isolate the tooth with a rubber dam.
- Remove the infected dentine with a round bur and/or a hand excavator. Leave the affected dentine.
- Adapt a matrix around the tooth if a wall is missing.
- Prepare Biodentine™ as indicated above (Biodentine™ mixing instructions)
- Insert Biodentine™ in the cavity avoiding trapping air bubbles. Ensure good adaptation of the material to the cavity walls and margins. Do not apply excessive pressure on the material.
- Model the surface of the restoration
- Wait until the end of the setting time before removing the matrix
- To optimize the mechanical properties of the material and facilitate removal of the matrix, a varnish can be applied onto the surface of the restoration.
- Check occlusion.
- Within one week to six months after placement of Biodentine™, prepare the cavity according to the criteria recommended for the selected restorative material.

The remaining Biodentine™ material can be considered as sound artificial dentine and permanently left in deep areas of the cavity and in areas adjacent to the pulp chamber. Biodentine™ is compatible with all direct or indirect crown restoration techniques (Inlay/Onlay), and particularly with all types of bonding systems.

##### PULP CAPPING (DIRECT AND INDIRECT):

Assess pulp vitality by the usual tests.

- Isolate the tooth with a rubber dam.
- Remove the infected dentine with a round bur and/or a hand excavator. Leave the affected dentine.
- Adapt a matrix around the tooth if a wall is missing.
- If there is bleeding in the pulp, hemostasis must be achieved before applying Biodentine™.
- Prepare Biodentine™ as indicated above (Biodentine™ mixing instructions).
- Place Biodentine™ directly on the exposed pulp avoiding trapping air bubbles. Ensure good adaptation of the material to the cavity walls and margins. Do not apply excessive pressure on the material.
- Perform the immediate or non-immediate enamel restoration as indicated above.
- Patients should be followed according to current recommendations.

##### PULPOTOMY:

Assess pulp vitality by the usual tests. In case of clinical signs and symptoms of irreversible pulpitis, pulpotomy is recommended when bleeding can be controlled in 5 minutes.

- Isolate the tooth with a rubber dam.
- Remove the infected dentine with a round bur and/or a hand excavator.
- Gain access to the pulp chamber and clean out the pulp.
- If there is bleeding in the pulp, hemostasis must be achieved before applying Biodentine™. If haemostasis cannot be achieved after 5 minutes, further pulp tissue should be removed (partial or full pulpotomy) step by step until a controlled bleeding. A full coronal pulpotomy can be carried out to the level of the root canal orifices with bleeding arrested.
- Adapt a matrix around the tooth if a wall is missing.
- Prepare Biodentine™ as indicated above (Biodentine™ mixing instructions).
- Place Biodentine™ directly in the pulp chamber and ensure good adaptation to the cavity walls and margins
- Model the surface of the restoration.
- Wait until the end of the setting time of the material before removing the matrix.
- To optimize the mechanical properties of the material and facilitate removal of the matrix, a varnish can be applied onto the surface of the restoration.
- Check occlusion.

- Within one week to six months after placement of Biodentine™, prepare the cavity according to the criteria recommended for the selected restorative material.
- Patients should be followed according to current recommendations. The remaining Biodentine™ material can be considered as sound artificial dentine and permanently left in deep areas of the cavity and in areas adjacent to the pulp chamber. Biodentine™ is compatible with all direct or indirect crown restoration techniques, and particularly with all types of bonding systems.

##### REPAIR OF ROOT PERFORATIONS:

- Isolate the tooth with a rubber dam.
- Prepare the root canal alternately using suitable endodontic instruments and a solution of sodium hypochlorite.
- Dry the canal with paper points without totally dehydrating the root canal and use a chlorhexidine solution or a calcium hydroxide paste for disinfection between visits. Tightly seal the access cavity with a temporary cement to protect the temporary filling.
- At the next visit (usually after one week), place a rubber dam and remove the temporary crown restoration. Clean the canal alternately using a solution of sodium hypochlorite and suitable endodontic instruments. Dry the canal with paper points without totally dehydrating the root canal.
- Prepare Biodentine™ as indicated above (Biodentine™ mixing instructions).
- Dispense Biodentine™ over the perforation site using a suitable instrument.
- Condense Biodentine™ with a plugger.
- Take an X-ray to check that the material is correctly positioned.
- Remove excess material and place a temporary filling.
- Complete root canal treatment at the next visit according to current recommendations.

##### REPAIR OF FURCATION PERFORATIONS:

- Isolate the tooth with a rubber dam.
- Rinse the cavity with a solution of sodium hypochlorite to disinfect the area.
- If there is bleeding, hemostasis must be achieved before applying Biodentine™.
- Dry the pulp chamber.
- Prepare Biodentine™ as indicated above (Biodentine™ mixing instructions).
- Dispense Biodentine™ and condense. Perforation repair and crown restoration are performed in a single step.
- Take an X-ray to check that the material is correctly positioned.
- Remove excess material.
- At a subsequent visit, if all clinical signs of a successful treatment are present, the possibility of a permanent restoration can be considered.

##### REPAIR OF PERFORATING INTERNAL RESORPTIONS:

- Isolate the tooth with a rubber dam.
- Prepare the root canal alternately using suitable endodontic instruments and a solution of sodium hypochlorite.
- Dry the canal with paper points without totally dehydrating the root canal and use a calcium hydroxide paste for disinfection between visits. Tightly seal the access cavity with a temporary cement to protect the temporary filling.
- At the next visit (usually after one week), place a rubber dam and remove the temporary crown restoration. Clean the canal alternately using a solution of sodium hypochlorite and suitable endodontic instruments. Dry the canal with paper points without totally dehydrating the root canal.
- Prepare Biodentine™ as indicated above (Biodentine™ mixing instructions).
- Dispense Biodentine™ over the resorptive defect using a suitable instrument.
- Condense Biodentine™ with a plugger.
- Take an X-ray to check that the material is correctly positioned.
- Remove excess material and place a temporary filling.
- Complete root canal treatment at the next visit according to current recommendations.

##### APEXIFICATION:

- Isolate the tooth with a rubber dam.
- Prepare the root canal alternately using suitable endodontic instruments and a solution of sodium hypochlorite.
- Dry the canal with paper points without totally dehydrating the root canal and use a calcium hydroxide paste for disinfection between visits. Tightly seal the access cavity with a temporary cement to protect the temporary filling.
- At the next visit (usually after one week), place a rubber dam and remove the temporary crown restoration. Clean the canal alternately using a solution of sodium hypochlorite and suitable endodontic instruments. Dry the canal with paper points without totally dehydrating the root canal.
- Prepare Biodentine™ as indicated above (Biodentine™ mixing instructions).
- Dispense Biodentine™ in the root canal using a suitable instrument.
- Condense Biodentine™ with a plugger.
- Take an X-ray to check that the material is correctly positioned.
- Remove excess material and place a temporary filling.
- Complete root canal treatment at the next visit according to current recommendations.

##### ROOT-END FILLING IN ENDODONTIC SURGERY:

- Gain access to the operative site following the current recommendations in endodontic surgery.
- Using a specific ultrasonic tip, prepare a root-end cavity, 3 to 5 mm deep in the apical portion of the root canal.

- Isolate the area. Achieve hemostasis. Dry the cavity with paper points.
- Prepare Biodentine™ as indicated above (Biodentine™ mixing instructions).
- Dispense Biodentine™ in the cavity using a suitable instrument. Condense Biodentine™ with a small plugger.
- Remove excess material and clean the surface of the root.
- Take an X-ray to check that the material is correctly positioned.

##### REVITALIZATION

###### First appointment

Assess pulp vitality and apex diameter by the usual tests. Revitalization procedure is indicated for an immature tooth with necrotic pulp.

- Isolate the tooth with a rubber dam.
- Remove loose or necrotic pulp tissue using suitable endodontic instruments. Avoid mechanical instrumentation of the root canal walls.
- Irrigate with sodium hypochlorite solution (advised concentration: 1.3 to 3 %, 20 mL, 5 min), use of side-vented needle, place 2 mm above vital tissue.
- Irrigate with sterile physiological saline (5 mL) and dry with paper points. Dry then with EDTA (advised concentration: 15 – 17 %, 20 mL).
- Insert a non-discoloring calcium hydroxide product homogenously into the root canal. Instead of calcium hydroxide, a triple antibiotic paste consisting of ciprofloxacin, metronidazole and minocycline can be used.
- Place coronal seal directly onto canal dressing with a minimum thickness.

###### Second appointment (usually 2 to 4 weeks later)



- Anesthesia with a local anesthetic without vasoconstrictor. Remove the temporary seal.
- Irrigate with EDTA (advised concentration: 15 – 17 %), use a side-vented needle, place 2 mm above vital tissue. Irrigate with sterile physiological saline (5 ml). Remove excess liquid with paper points.
- Induce mechanical bleeding of periapical tissue and rotational movement of an apically pre-bent file. Allow the canal to fill blood until 2 mm below the gingival margin, wait for blood clot formation for approximately 15 minutes.
- Prepare Biodentine™ as indicated above (Biodentine™ mixing instructions)
- Cut a collagen matrix to a diameter larger than the coronal part of the root canal and a height of 2 - 3 mm, place on top of the blood clot, allow the matrix to soak with liquid, avoid formation of a hollow space.
- Place Biodentine™ on top of the collagen matrix in a thin homogenous layer of about 2 mm underneath the cement-enamel junction.
- Then fill the cavity with Biodentine™.
- On a next visit, refresh the cavity walls with a diamond bur or grit blast with aluminum oxide and seal with adhesive restoration.
- Patients should be followed according to current recommendations (6, 12, 18 and 24 months).

#### STORAGE AND HANDLING CONDITIONS

Do not store above 25°C and protect from moisture.

Use the 5-liquid single-dose containers within 3 months after first opening of the foil.

Liquid in single-dose containers must be used immediately after opening. Any remaining material shall be discarded.

	Upper limit of temperature
	Consult instructions for use
	Keep dry, protect from excessive humidity, rain or damp conditions
	Do not reuse
	Catalogue or model number
	Batch code
	Use-by date
	Manufactured by