




# COMPACT *SLEEK* DISCREET



## PAROJECT® INTRALIGAMENTAL SYRINGE

Paroject is an intraligamental syringe used for “PDL” injections. This pen style syringe is compact, durable and most of all - easy to use. Its unique design makes it versatile and less intimidating to the patient.

- Permanent hub design
- Unintimidating in appearance
- Can be sterilized by any method
- Made from high-grade stainless steel
- Injects only 0.06 ml. per “click”, permitting a slow non-traumatic injection
- Uses standard disposable dental needles
- Compact and user-friendly (delicate pen-type design)



### Description

Injection of local anesthetic solution

### Product

### Item #

Aspiject Intraligamental Syringe

01N2040


To order Septodont Paroject syringe, contact your dental dealer or call Septodont at 800-872-8305

800-872-8305 • septodontusa.com  
From the manufacturers of Septocaine®  
anesthetics • materials • endodontics • infection control





# BALANCE *FEEL* CONTROL



## ASPIJECT® SELF-ASPIRATING SYRINGE

Injection of local anesthetic solution

- Can be sterilized by any method
- Made from high-grade stainless steel
- Uses standard disposable dental needles
- Well balanced and comfortable use
- Lightweight and durable while providing excellent tactile sensitivity during local anesthesia injections



### Description

Injection syringe for infiltration and block anesthesia

### Product

### Item #

Aspject Self-Aspirating Syringe - Saddle grip  
 Aspject Self-Aspirating Syringe - Thumb grip

01N2020  
 01N2030

To order Septodont Aspject syringes, contact your dental dealer or call Septodont at 800-872-8305

800-872-8305 • septodontusa.com  
From the manufacturers of Septocaine®  
anesthetics • materials • endodontics • infection control

